# INSTRUCTION MANUAL

## Digital and Optical Microscope

Model # 44330


## Introduction

Congratulations on your purchase of a Celestron microscope. Your microscope is a precision optical instrument made of high quality materials to ensure durability and long life. It is designed to give you a lifetime of pleasure with a minimal amount of maintenance.

Before attempting to use your microscope, please read through the instructions to familiarize yourself with the functions and operations to make usage easier for you. See the microscope diagram to locate the parts discussed in this manual.

The microscope provides high powers from 100x to 400x using it as an optical microscope and much higher powers when used as a digital microscope. It is ideally suited for examining specimen slides of yeasts and molds, cultures, plant and animal parts, fibers, bacteria, etc. You can also examine small and thin objects at the lower powers such as coins, rocks, insects, circuit boards, various materials, and many other objects.

With the included digital camera and the software, you can observe the magnified images or capture video or take snapshots.

**Figure 1**

## Standard Accessories with your Microscope


|                      |                               |
|----------------------|-------------------------------|
| • 20x Eyepiece       | • Digital Camera w/ USB Cable |
| • 5x Objective Lens  | • 3 Prepared Slides           |
| • 10x Objective Lens | • CD-ROM – ArcSoft Software   |
| • 20x Objective Lens | • Instruction Manual          |
| • Bottom Illuminator |                               |

## Specifications

| Model # 44330                          | Specifications                                                      |
|----------------------------------------|---------------------------------------------------------------------|
| Stage                                  | Plain Stage with Spring Loaded Clips – 80 mm x 60 mm (3.1” x 2.4”)  |
| Eyepiece                               | Glass Optics                                                        |
| Focuser                                | Coarse Focus – Dual Knobs                                           |
| Objective Lenses                       | All Glass Optics – 5x (100 power), 10x (200 power), 20x (400 power) |
| USB Cable                              | 2.0                                                                 |
| Illuminator -- Bottom                  | LED -- Uses 2AAA Batteries (user supplied)                          |
| Nosepiece                              | Triple with Click Stop                                              |
| Digital Camera                         | 2MP CMOS – 1600x1200 pixel array                                    |
| Weight (with batteries) and Dimensions | 19oz (539g) 4.5” x 4” x 6” (114 mm x 102 mm x 152 mm)               |

## Setting Up Your Microscope

1. Carefully remove the microscope and other parts from the box and set them on a table, desk, or other flat surface.
2. Remove the plastic bag covering the microscope.
3. Remove the plastic plug up cap from the eyepiece tube and install the eyepiece – see Figure 1.
4. Install the batteries for the bottom illuminator in the base of the microscope which uses two AAA batteries (user supplied). See Figure 2 showing the open battery compartment under (bottom) microscope. Install the batteries (see the inside of the battery compartment where the positive (+) and negative (-) ends of the batteries go).


**Figure 2**

**You are now ready to use your microscope for looking at specimen slides or small objects through the microscope eyepiece as an optical microscope! For using the digital (CMOS) camera, installation and operating instructions will follow later in this manual.**

# Microscope Operation

## Viewing a Specimen Slide

**Figure 3**

Carefully place a specimen slide under the stage clips and center the specimen directly over the hole in the center of the stage (where the bottom illuminator light shines up through a specimen slide – see Figure 3 to the left. It will take some experimenting to place slides or objects in the center of the stage as the image you see is upside down and reversed but after some usage you will have an easy time centering. Read the sections below on Focusing, Changing Power, and Illumination before proceeding.

**Viewing Objects** – you can view coins, stamps, currency, plants, foods, circuit boards, and many other things. You should use only use 100x and 200x powers and the height should be no more than ¼” (6 mm). You will also need to use a bright desk lamp (or similar) that can focus directly on the object. A well slide would be helpful to center very small objects.

## Focusing & Changing Power (Magnification)

Now that the specimen slide (or object) is placed directly under the objective lens, use the focus knob to focus on the specimen. The focus knob takes some getting used to – rotate past the focus and then back as you will see how to obtain the best focus position. After some usage the sharp focus position will be easy to find.


1. Always start with the lowest power (5x objective lens) which will give you 100 power.
2. For higher powers, you will have to rotate the nosepiece to change the objective lens to the 10x one (200 power) or the 20x one (provides power of 400x). You rotate the nosepiece by holding the microscope above the nosepiece with one hand and rotate the nosepiece with the other hand until it clicks at the position.  
**Note: The power of the objective lens you are using is shown on the nosepiece after it clicks into position.**
3. At the highest powers, your views will be greatly magnified but somewhat darker. The most enjoyable views can be at the lower powers which have a wider field of view and brighter illumination.

**Note: You can obtain other powers if using optional 23.2 mm optional eyepieces obtained from microscope accessory supply companies.**

## **Illumination**

To get the sharpest and best views, the illumination (lighting) will have to be adjusted.

1. The bottom illuminator is used for specimen slides which shine up through the hole in the stage through the slide. The illuminator is turned on by a switch dial shown in Figure 4.
2. The light from the bottom illuminator can be increased or decreased by moving the illuminator dial. You will have to experiment to provide the best lighting for the best views.
3. To get the illumination necessary for objects you will have to use a bright desk (or similar) lamp that focuses the light directly on the object. You will have to experiment to get the best light for your needs.

**Figure 4**

## **Using the Digital Camera for Viewing and/or Imaging with your Microscope**

### **Attaching your Digital Camera to your Microscope**

First, you need to remove the eyepiece (see Figure 1) from the eyepiece tube by pulling it straight up. Take off the protective cap from the camera. Next, install the camera into the eyepiece tube by pushing it gently into the eyepiece tube. Lastly, you plug the USB cable into an open USB port on your computer.

When viewing or imaging a specimen slide or object, you can change the orientation of the image on the computer screen by rotating the camera to the position you desire.

Left to right on next page images – digital camera with USB cable, eyepiece tube, and digital camera attached to eyepiece tube.

**Figure 5a**

**Figure 5b**

**Figure 5c**

Before using your microscope for viewing on your computer screen or imaging, you will have to ensure the driver is recognized by your computer. Your microscope is a “Plug and Play” type where the computer should automatically recognize the driver type on Windows XP and Vista (32 and 64) Operating Systems as well as MAC computers operating on Version 10.4.9 and higher.

Plug the USB cable of the digital camera into the USB port of your computer. Each computer system will have slightly different messages when you attach the digital camera. For example, on some Vista OSs when plugged in the message may say “Installing Device Driver Software” and shortly thereafter you may see a message “Your Devices Are Ready to Use”. Just follow any instructions shown on the screen.

You will need to install the ArcSoft software to actually view and/or image and prior to doing this you will need to unplug the USB cable. Note that if using a MAC computer you will need to have an image capture program for your MAC as ArcSoft will not work on MAC computers.

### **Installing the Viewing and Imaging Software**

Your microscope comes with the ArcSoft Software CD-ROM which allows you to use your microscope for visual purposes and for imaging with your computer.

Before installing the software, remove the USB cable of the Digital Camera from your computer and do not connect it until the software installation is complete. Insert the CD-ROM into your CD-ROM drive and it should automatically begin. Messages will appear (each operating system may show something somewhat different) to guide you through the installation process.

If the software did not begin automatically, double click on (My) Computer and right click the CD-ROM icon and it should start.

- You will see “Files Currently on the CD” and notice “ArcSoft” as one of the folders – other folders will be Web Companion, Help, 44306 Instructions, and 44330 instructions. Double click on “ArcSoft” and then among the choices double click on “Setup” (InstallShield Setup Launcher).

- You will be asked to choose which language you want to use. English is the default language with 13 other language choices.
- Then, the InstallShield Wizard will appear. Enter the License Key Number if asked for. Then select the default choices, Next, etc. until “Finish” is seen and click on it.

Your software was properly installed. You will have a desktop icon called “WebCam Companion” but if not that will be the program you select from your installed programs when you want to use the digital microscope.

The software allows you to observe specimen slides or objects on your computer. When you view with the camera installed, the magnification depends on the objective lens you are using and also the size of your PC monitor.

### **Magnification Estimated Using the Digital Camera**

|                                       |       |                         |       |
|---------------------------------------|-------|-------------------------|-------|
| 5x Objective using a 14” monitor ---  | 240x  | using a 17” monitor --- | 290x  |
| 10x Objective using a 14” monitor --- | 850x  | using a 17” monitor --- | 1030x |
| 20x Objective using a 14” monitor --- | 1550x | using a 17” monitor --- | 1880x |

### **Using the ArcSoft Software**

The software contains many applications for capturing and processing your photos (snapshots) or videos. Click on WebCam Companion to begin.

The Application Panel is very intuitive. The screen will be divided into “My Captures” and “Capture”. The screens have various buttons and just put the cursor over them to see what they are and what they do. To begin you should go to “Preferences” and choose a selection in the “Media Folder” and then you can use the default or change settings for “Capture Settings” and “Web Camera Settings”.

If any questions or problems come up, you can use the “Help” menu.

Plug in the USB cable of the Digital Camera to your computer.

You will use the illuminator, focuser, and stage the same way as you would with the optical microscope operation.

***You are now ready to begin your digital adventure!***

## Care, Maintenance, and Warranty

Your microscope and digital camera are precision optical instruments and should be treated with care at all times. Follow these care and maintenance suggestions and your microscope will need very little maintenance throughout its lifetime.

- When you are done using your microscope, remove any specimens or objects left on the stage.
- Turn off the illuminator when you are done using the microscope.
- If you will not be using your microscope for a long period of time, remove the batteries from the illuminator.
- Always place the dust cap on the eyepiece or on the camera sensor when not in use or when being stored.
- Store the microscope in a dry and clean place.
- Be very careful if using your microscope in direct sun light to prevent damage to the microscope or your eyes.
- Never point the sensor on the camera towards the sun or the camera can be damaged and cease working.
- When moving your microscope, carry it by the “arm” with one hand.
- Clean the outside surfaces with a moist cloth.
- Never clean optical surfaces with cloth or paper towels as they can scratch optical surfaces easily.
- Blow off dust with a camel’s hair brush or air blower off optical surfaces.
- To clean fingerprints off of optical surfaces, use a lens cleaning agent and lens tissue available at most photo outlets and when cleaning do not rub in circles as this may cause streaks or scratches to occur.
- Never disassemble or clean internal optical surfaces. This should be done by qualified technicians at the factory or other authorized repair facilities.
- Do not disassemble or take apart the digital camera or damage can occur.
- When handling glass specimen slides, use care as the edges can be sharp.

### Warranty

Your microscope has a two year limited warranty. Please see the Celestron website for detailed information on all Celestron microscopes at [www.celestron.com](http://www.celestron.com).

#### FCC Statement

This device complies with Part 15 of FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference, and
2. This device must accept any interference received, including interference that may cause undesired operation.

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Torrance, California 90503 U.S.A.