UNIVERSAL DIGITAL LCD WITH CAMERA FOR MICROSCOPES

INSTRUCTION MANUAL

MODEL #44430

ENGLISH

INTRODUCTION

Thank you for purchasing the Celestron universal LCD Monitor with Digital Camera Accessory (LMDC). This product is a precision optical instrument, made of the highest quality materials to ensure durability and long life.

Before attempting to use the LCD Monitor with Digital Camera, please read through the instructions to familiarize yourself with the functions and operations to maximize your enjoyment and usage.

The LCD Monitor and Camera does not use eyepieces. You will view specimens or objects on the LCD screen which are easy to see and you can enjoy them with others also. You can take

snapshots or short videos with the supplied digital camera. You can also view your images on a TV monitor.

The final section provides simple care and maintenance tips for you to follow to ensure that your accessory provides you with years of quality performance, usage, and enjoyment.

Note: The LCD monitor viewing will be full screen size, 3.5" (89 mm), when using an image tube of a microscope (if available) or a 2" (50 mm) circular image when using an eyepiece tube used for normal viewing with eyepieces.

STANDARD ITEMS AND SPECIFICATIONS

| | |
|----------------------------|---|
| LCD Monitor – 3.5" (89 mm) | 2" circular image when using eyepiece tubes and full screen viewing when using image tubes. High Definition Digital TFT display – 320x240 pixels |
| Digital Camera | 3MP CMOS; 10x magnification in lieu of an eyepiece; 2048x1536 pixel array; JPEG for snapshots; AVI files for video – 640x380 pixels (30fps) |
| Flash Memory – Internal | 32MB |
| SD Card Slot | Optional SD Card can be used up to 2GB |
| LCD Positioning | 100° vertical range (from even position – 60° up and 40° down; 360° horizontal – positioning on most microscopes) |
| Adapters/Tubes | 23.2 mm and 30 mm Tube Adapters; 12.5 mm Joint Adapter/Connector (CS type adapter); 17.5 mm Joint Adapter/Connector (C type adapter) |
| Projection Mirror Tube | |
| USB Cable | 1.1 Mini – length 5' (1.5m) |
| Instruction Manual | |
| Power/Current | Output – 5-Volt DC / 1A |
| AC Adapter | Input – Universal 100 to 240 Volt, 50/60HZ; Four plugs included – USA/EU/UK/Australia |
| Weight/Dimensions | 10oz (283g) / 5" (127 mm) x 6.5" (165 mm) x 3.5" (89 mm) |

SETTING UP YOUR LCD MONITOR AND DIGITAL CAMERA

1. Remove the case (Figure 1) containing the various parts from the carton.

2. Carefully remove the LCD monitor, Digital Camera and other parts from the case (Figure 2) and set them on a table, desk, or other flat surface.

There are additional parts which will be discussed later.


FIGURE 1

FIGURE 2

FIGURE 3


The LCD Monitor functions will be discussed in detail in the various setup sections. In Figure 3 is shown the various function keys for the LCD Monitor.

ATTACHING THE LCD MONITOR AND DIGITAL CAMERA TO VARIOUS TYPES OF MICROSCOPES

The LCD Monitor and Digital Camera are designed to work with virtually all traditional size microscopes (not toy microscopes) – those with 23.2 mm eyepiece tubes (biological/compound microscopes) and those with 30 mm eyepiece tubes (stereoscopic microscopes). It doesn't matter if the microscopes are monocular or binocular types. In addition, the LCD Monitor and Digital

Camera will adapt to most microscopes that have a 2nd or 3rd tube called by various names (teaching tube/image tube/extra eyepiece tube) and you will use one of the standard type threaded joint adapters/connectors of 12.5 mm (CS type) or 17.5 mm (C type) sizes along with various adapters from the microscope manufacturer.

FIGURE 4 – ADAPTATIONS TO VARIOUS MICROSCOPE TYPES

ATTACHING THE LCD MONITOR AND DIGITAL CAMERA TO STANDARD EYEPiece TUBES – SEE THE FIRST THREE TYPES IN FIGURE 4

1. In Figure 5a is shown the back of the LCD Monitor and under the round part is the Digital Camera. Remove the protective cap from the round part which protects the digital camera.
 2. In Figure 5b is shown where the Projection Mirror Tube will thread into the Digital Camera. Remove the protective caps from both ends of the projection mirror tube. Then, thread the projection mirror tube (the large end of the tube) clockwise into the threads of the digital camera. Figure 5c shows the projection mirror tube attached to the digital camera.
 3. Remove the standard eyepiece from the eyepiece tube (if using a monocular microscope there is only one eyepiece to remove). If using a binocular microscope (using two eyepieces), you can remove the eyepiece from either eyepiece tube.
 4. For biological microscopes, you will use the 23.2mm Tube Adapter (the 30mm Tube Adapter would be used on stereo microscopes). Figure 5d shows the 23.2mm Tube Adapter on top and it will be threaded clockwise onto the smaller end of the projection mirror tube until tight.
 5. The LCD Monitor with Digital Camera and Projection Mirror Tube attached can now be inserted into the eyepiece tube of the microscope and it is held in place by a friction fit.
- Note: Be careful when moving your microscope with the LCD Monitor and Digital Camera attached. Keep the LCD Monitor upright so that it will not accidentally fall out.*

FIGURE 5A

FIGURE 5B

FIGURE 5C

FIGURE 5D

ATTACHING THE LCD MONITOR AND DIGITAL CAMERA TO IMAGE TUBES

Image Tube adaptations are quite varied as different microscope manufacturers offer different types of adapters. One such configuration is shown in the far right of Figure 4.

Typically you will use one of the threaded joint adapters which will thread into the digital camera – Figure 5a. Then, you will use

adapters from the particular microscope manufacturer that will be made for your microscope. Some manufacturers use simple tube adapters with built-in lenses and others use complex and several different tubes, etc.

OPERATION OF THE LCD MONITOR AND DIGITAL CAMERA WITH YOUR MICROSCOPE

Before looking at slide specimens or objects you must remove the protective cover over the LCD screen, turn the LCD on, turn on the proper illuminator of your microscope, focus as you normally do, and understand how to use the LCD Monitor with Digital Camera.

All functions of your microscope are operated the same as you always have done with the only difference being that you will be using the LCD monitor with its attached digital camera for viewing and imaging. The digital camera replaces the eyepiece(s) with a 10x magnification itself and the power of the microscope is then determined by the objective lens chosen times 10.

FIGURE 6

FIGURE 7

TO POWER UP THE LCD MONITOR

1. Plug the small cable from the AC adapter into the socket on the side of the LCD Monitor (bottom of Figure 6) labeled 12V DC.
2. Insert the plug end of the AC adapter into the proper power source.

Note: The AC Adapter supplied with this microscope has an interchangeable plug system that can work in the USA, UK, Europe, Australia, and many other parts of the world. The plug for the USA will be installed in the AC Adapter. You can easily change plugs for your particular requirement by doing the following:

- A. Pull the plug up and away while holding the base of the AC adapter and the plug will come off.

- B. Take the plug style that you want to use and align the two holes in the bottom of the plug with the two pins in the AC adapter (Figure 7). Then, push the plug into the AC adapter all the way and you're finished.

Prior to viewing specimens or objects with your microscope, you need to turn the LCD Monitor on by pushing the Power On key (Figure 3). Push the button and shortly the Monitor On LED (green color) will come on indicating the LCD Monitor is "On". That is all you need to do to use the LCD screen for viewing. The various keys on the LCD Monitor are mainly for taking images (snapshots and video) and will be discussed later in this manual.

DIGITAL IMAGING

You can take snapshots or a short video with your microscope with the built-in digital camera. With the internal storage memory, you do not need to use a PC or any other devices to do imaging. Transferring of images to your PC for saving and or printing them is easy and will be discussed later in this manual.

Note: If you are going to take images, do not connect the USB cable to your PC or damage could occur to the equipment. The USB cable is not used at all for taking images.

SETTINGS AND INFORMATION FOR THE DIGITAL CAMERA

1. In the bottom right of the LCD screen is the EV (Exposure Value) which adjusts the brightness level. The normal position is 0.0 and this can be adjusted in increments up to +1.2 to -1.2 by pushing the up/down directional keys.
2. In the bottom left of the screen is shown an estimate of the remaining snapshot images left in the internal memory.
3. In the top left of the screen is shown an icon of a camera which indicates it is ready for snapshot mode.
4. By pushing the directional keys (left/right) you change the 4x digital zoom shown by a magnifier icon when you push the keys. The digital zoom is enabled when using resolution settings of 1024x768, 800x600, or 640x480 and does not function at higher resolution settings.
5. Push the Menu Key to:
 - A. Mode – the normal (default setting) is "Single" for taking snapshots one at a time as you choose. You can change this setting to "Continuous" which will tell the camera to take snapshots at a predetermined time you set in the "Setting" choice.
 - B. Size – this sets the resolution you want for the image. The camera itself uses the 2048x1536 pixel arrangement (3MP) but you can use digital interpolation for lower

resolutions to increase the number of images you can save in your internal storage memory. Use the directional keys (up/down) to make changes and you will see the resolution setting at the top left of the screen.

- C. Effect – you can change from Normal, to Black/White, Negative, or Sepia.
- D. Date Label – if you desire to have the time and or date stamped on your images you can select the information to include under "Setting" and then choose "Yes", or you can choose to have no label (default setting).
- E. Video – if you want to take videos, push the OK key and the resolution will appear on the top left of the screen. Resolution is automatically changed to low resolution for video. The rate is 25 fps. The time elapsed for the video is shown in the lower left of the screen. To begin a video push the "SNAP" button and end the video by pushing the "SNAP" button. While recording a video icon will also be flashing in the top left of the screen. To exit to the snapshot mode, push the "ESC" button.

Note: After turning the LCD Monitor off, most settings will revert to the default setting.

TAKING IMAGES




1. Snapshots – to take snapshots make sure you are in the camera mode (camera icon in the lower left of the screen) and you will use the Snap Key on the LCD Monitor. Push the key when ready and the image will be stored in the internal storage memory.
2. Video – change the settings to the Video. To begin and end your video recording push the "SNAP" key. The movie (video) icon will flash to let you know you are recording and the elapsed time will be shown in the bottom left of the screen.

3. You can also take images (snapshot or video) by using an SD (Secure Digital) Card. Your SD card is inserted in the LCD Monitor (1). When you insert the SD card properly the SD card icon will be visible at the top of the LCD screen. On the bottom right of the screen will be shown the memory available on the SD card. When imaging this way the images are saved on the SD card and not in the internal storage memory.

Note: Inserting or removing an SD card while the LCD is on may cause the LCD to shut down and/or could damage the SD card.

SAMPLE IMAGES

The following specimen images left to right below were taken at 40x, 100x, 400x of Vicia Faba Young Root C.S.

REVIEWING AND MANAGING YOUR IMAGES

You can review and manage your snapshots and videos on the LCD Monitor that are saved in the internal storage memory.

1. By pushing the ESC Key, you enter the stored images (Picture View). Use the directional keys to select the image you want to view and press the OK Key to view with full screen. Thumbnail images appear on the left half of the screen (as you scroll through the images on the right) and the images are listed by number on the right side of the screen for ease of selection. To review your videos while in the stored images, press MENU and then select the video you want to review. Push OK to start and stop the video. Press ESC to exit the video. When observing an image with full screen view, you can return to

the stored images by pushing the ESC Key. To exit to real time (snapshot mode), push the ESC Key again.

2. You can manage your images by pushing the MENU Key:
 - A. File Protect – you can make selections by pushing the OK Key. You have a choice of "Lock" to protect this image or "Unlock" where the image is not protected.
 - B. Del File – if you want to delete one or more images you can select "Current" to delete the image you are viewing or you can select "All" to delete all of your images.
 - C. Exit – select this to return to the "Picture View" screen.

TRANSFERRING YOUR IMAGES

To transfer images to a PC or MAC, you need to have a free USB port and have an imaging program (for snapshots and/ or video). Turn the power Off of the LCD Screen.

Note: Do not disconnect the USB cable while transferring images or damage may occur.

1. You can transfer images from the internal storage memory to your PC by using the supplied Mini USB Cable (located in a pullout drawer of the case). The small plug end of the cable plugs into the LCD Monitor and the large plug end of the cable plugs into your PC. Turn the power On of the LCD screen. If the connections are proper you will see on

the LCD screen "MSDC" or similar data. Your PC will automatically recognize the new hardware. Then, you will choose which program on your PC you want to transfer the images to – you may have to choose a different snapshot program and a different video program depending on what software imaging programs you have.

2. If you used an SD card to store your images, you can transfer them to your PC using the method in # 1 above or you can take the SD card out and use the SD card slot on your PC (if available) to transfer the images. In either way your PC will ask you to choose which program you want to transfer the images to.

VIEWING IMAGES ON A TV SCREEN

Video cables are supplied (in the drawer of the case). Plug one end of the video cable (red and white connectors) into the TV set "Video In" and "Audio In" sockets. Plug the other end of the cable (black) into a socket on the side of the LCD Monitor labeled "Video Out". Press MENU and then select "TV Output" and

choose LCD Output, NTSE Output, or PAL Output. The image as seen on your LCD screen (in the Preview mode) will now be viewable on your TV screen. The image on the TV screen will be greatly enlarged and not as sharp in color or contrast as the LCD screen due to the resolution of the TV screen.

CARE, MAINTENANCE, AND WARRANTY

Your Celestron accessory is a precision optical instrument and should be treated with care at all times. Follow these care and

maintenance suggestions and your microscope will need very little maintenance throughout its lifetime.

- When you are done using your microscope, remove any specimens left on the stage.
- Turn off the illuminator switch.
- Turn off the LCD monitor.
- Unplug the power cord(s).
- Remove the LCD Monitor and Digital Camera from your microscope and place them back in the case.
- Store the product in a dry and clean place.
- When moving your microscope, make sure the LCD Monitor is in the upright position and cannot accidentally fall.
- Clean the outside surfaces with a moist cloth.
- Always unplug any cords before cleaning.
- Never clean optical surfaces with cloth or paper towels as they can scratch optical surfaces easily.
- Blow off dust with a camel's hair brush or an air blower from optical surfaces.
- To clean fingerprints off of optical surfaces, use a lens cleaning agent and lens tissue available at most photo outlets and when cleaning do not rub in circles as this may cause sleeks or scratches to occur.
- Never disassemble or clean internal optical surfaces. This should be done by qualified technicians at the factory or other authorized repair facilities.
- When handling glass specimen slides, use care as the edges can be sharp.

**YOUR CELESTRON ACCESSORY HAS A TWO YEAR LIMITED WARRANTY.
PLEASE SEE THE CELESTRON WEBSITE FOR DETAILED INFORMATION ON ALL
CELESTRON ACCESSORIES AT WWW.CELESTRON.COM.**

EEC: This product complies with EEC guidelines in EN61558-2-6:1997 and EN61558-1:1997+A1

FCC Statement

This device complies with Part 15 of FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference, and
2. This device must accept any interference received, including interference that may cause undesired operation.

2835 Columbia Street • Torrance, CA 90503 U.S.A.

Telephone: 310.328.9560 • Fax: 310.212.5835

©2010 Celestron

All rights reserved. • Printed in China • 09-10

Product design and specifications are subject to change without prior notification.

Designed and intended for those 13 years of age and older.

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